<b>VILLACRYL HARD</b>	<b>INSTRUCTIONS FOR USE - EN</b>
-----------------------	----------------------------------

**I. Product description**

**Villacryl Hard** : Hard denture relining material.

**Processing parameters**: necessary parameters for the proper processing of the materials are given in further section of IFU.

**Indications for use**

**Villacryl Hard**: relining of acrylic removable complete and partial dentures directly in the patient’s mouth or indirectly in the dental laboratory based on the impression of the denture.

**Clinical benefits**

Restoring lost biting and chewing functions to users of dental prostheses that have lost their retention to the denture foundation.

**Target group or target groups of patients and intended users of the devices**

**Users**: Professionals in the dental technology industry.

Qualified dentists

Qualified dental technicians producing dental prostheses.

**Patients**: People of different ages and sexes, denture users that have lost their retention to the denture foundation.

**Contraindications**

Do not use the product by people diagnosed with a lack of resistance to allergenic factors, including the ingredients of the product. In such situations, use a different material suitable for the patient.

For people sensitive to taste or small sensations, with gag reflexes, it is recommended to perform indirect relining in the dental laboratory.

Do not use the product by pregnant women and nursing mothers in the group of dental technicians.

Do not reline in case of mucositis.

**Product form**

Separate: a mixture of powder in plastic containers, liquid and etchant in bottles.

The product is available in the form of a set: powder + liquid + etchant + hardener

**Chemical composition**

Powder: acrylic polymer, barbituric acid derivative, pigments, benzoyl peroxide,

Fluid: acrylic monomers, N,N dimethyl-p-toluidine (DMPT), hydroquinone
Etching agent: acetone, ethyl acetate

**Product compatibility**

For use with acrylic dentures.

**Safety instructions**

**Warnings**

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation

H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled

H335 May cause respiratory irritation.

H336 May cause drowsiness or dizziness

H411 Toxic to aquatic life with long-lasting effects

EUH066: Repeated exposure may cause skin dryness or cracking

The liquid is sensitive to temperature and UV radiation. An uncontrolled polymerization may occur.

**Precautions**

In case of any allergic reaction, stop work immediately and consult a dermatologist.

It is absolutely recommended that the patient rinse his mouth with water after relining, directly in the mouth.

In the case of patients who are sensitive to the components of the product, special care should be taken, do not reline directly in the patient’s mouth, and in special cases use a different material.

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P303 + P361 + P353 IF ON SKIN: Remove/Take off immediately all contaminated clothing and rinse skin with water/shower.

P304 + P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Contact an ophthalmologist is advisable in case of irritation.

P312 Call a POISON CENTER or doctor/physician if you feel unwell.

P403 + P233 Store in a well-ventilated place. Keep container tightly closed.

P233 Keep cool.

P501 Dispose of the contents/container to an authorized consignee of the waste.

It is recommended to use gloves made of HDPE polyethylene for hand protection, a filtering fluidcompatible face mask for respiratory protection, goggles and workwear.

Do not eat or drink while working with the material. In the event of accidental ingestion of an ingredient of the product, consult a doctor immediately.

Do not pour the liquid down the sanitation system.
Before using the product, check if the fluid retains its properties. In the event of visible traces of density change, cloudiness or discoloration, stop using the product.

**Information to be provided to patients / prosthesis users:**

To reduce the residual monomer, it is advisable to during the first three days, soak the relined prosthesis in water when not in use.

If any symptoms of an allergic reaction or irritation occur, discontinue using the prosthesis immediately and consult a doctor.

The color of the relined layer after a few months of use of the denture may change due to the consumption of food and drinks, pharmaceutical drugs or the use of chemical agents for disinfection and cleaning of the prostheses. It is a natural phenomenon; it does not pose a threat to the patient and does not reduce the usefulness of the prosthesis.

Strong oxidizing compounds contained in cleaning agents and mouthwashes may change the color of the prosthesis after polymerization. After each meal, the prosthesis should be washed with a soft brush under cool, running water.

After each meal denture should be cleaned with soft brush under cool tap water. For cleaning can be used unscented soap, various preparations for cleaning and caring for dental prostheses that do not contain abrasives. It is recommended to wash the dentures from time to time in a vibration washer filled with a special cleaning liquid for dentures. Report any damage of the denture to the dental technician. Repairs may be only performed by a dental technician. If a piece of damaged dentures is swallowed, consult a doctor immediately. The material after polymerization is not visible on X-ray pictures.

**Side effects**

Dental technicians/dentists

Irritation, redness or symptoms of hypersensitivity, including the occurrence of anaphylactic shock, in case of allergy to any of the product ingredients.

**Patients**

Direct relining technique:

A slight burning sensation in the mouth.

Irritation, inflammation and allergic reactions caused by the contact of the chemical components of the product with the mucosa.

In people sensitive to taste and small sensations, gag reflexes may occur.

**Assortment\***

**Color shades**

Villacryl Hard	
pink	

	<b>Powder + liquid + etchant + hardener kit</b>
Villacryl Hard	60 g powder + 40 ml liquid +10 ml etchant + 40 g hardener

**Problem prevention (1)**

See website: [www.everall7.pl](http://www.everall7.pl)

Please send any comments regarding our products to our address.

Information provided in any way, even during the demonstration, does not invalidate the instructions for use. Users are responsible for checking that the product is suitable for the intended use. The manufacturer is not liable for damages, including to third parties, resulting from non-compliance with the instructions or non-use. The manufacturer’s liability is in any case limited to the value of the products delivered.

**MD**

The symbol means a medical device as defined in EU Regulation 2017/745.

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